



BeneHeart R3 Vet

Veterinary Electrocardiograph

Your Faithful Consultant for Resting ECG Diagnosis

companion animal



Search Mindray Animal Medical



LinkedIn



Facebook



Instagram

www.mindrayanimal.com
service@mindrayanimal.com


©2026 Shenzhen Mindray Animal Medical Technology Co., LTD. All rights reserved.
P/N:ENG-BeneHeart R3 Vet-210285X4P-20260123



Reliable Results

- A 1,000Hz high digital sampling rate and 16,000Hz for pacemaker signal detection.
- Together with different filter setting, BeneHeart R3 Vet provides a reliable ECG waveform record comes with diagnostic results.

Unique Recorder




Compatible with both rolling paper and Z-fold paper, you can easily switch between these two styles of papers without dismantling the pressure lever.



Long Operation Time

Large capacity rechargeable Lithium ion battery is provided, enough for 500 records printing or a 6-hour continues work without printing.



Clear Display

5-inch color screen come with a high resolution of 800*480 pixels, offering a clear view of ECG waveforms with detail recorded.

Great Mobility

BeneHeart R3 Vet weighs only 1.2kg with battery, clinicians can easily carry the machine with single hand and perform resting ECG test in different scenarios.



Informatic Solution

Patient data can be saved as different formats, and transferred to PC or network server for historical data backup, review, printing and tele-diagnosis.

